



## APPLICATION BULLETIN #10

### DENTURES, CROWNS, BRIDGES

Newly-constructed dentures, crowns, and bridges require thorough cleaning before use. Flux, solder residue, buffing compounds, tripoli, rouge, plaster and stone, and other miscellaneous contamination should be cleaned completely. Cleaning with L&R Ultrasonics provide the fastest, most thorough cleaning possible.

Simply –

1. Fill the tank 3/4 full with General Purpose Cleaning Solution Non-ammoniated at the recommended dilution (see product label).
2. Place item into the Basket and cavitate for 5-10 minutes.
3. Rinse thoroughly with water.
4. Other solutions for specific contaminations are available. Follow directions on each container closely. Dentures that have tartar and other build-up should be cleaned using Tartar, Light Stain & Permanent Cement Remover. Place denture in Beaker with enough solution to cover it. Suspend Beaker in the ultrasonic machine using the Positioning Cover. The Beaker must be immersed at least one inch into the solution in the tank of the machine to ensure good acoustic coupling. Cavitate for 5-10 minutes and rinse denture with running water.

#### EQUIPMENT NEEDED:

Machines, Accessories\*: Q140 Ultrasonic Cleaning System with Timer, Positioning Cover, Beakers, Bur Tray, Auxiliary Pan, Basket

Solutions: General Purpose Cleaner Non-Ammoniated(PC228)  
Plaster & Stone Remover (PC230)  
Tartar, Light Stain & Permanent Cement Remover(PC232)  
Tartar, Light Stain Remover Powder (PC028)  
Alginate Remover (PC240)  
Advanced Formula Temporary Cement Remover (PC293)  
UltraDose Ultrasonic Cleaning Solution (PC012)  
UltraDose Enzyme Plus (PCUD038)

\*Recommended Ultrasonic Cleaning System varies as to internal tank dimensions and volume of required cleaning. Please contact an L&R Sales Specialist to discuss your specific cleaning needs.



*because clean matters*

L&R Manufacturing Company  
577 Elm Street, PO Box 607  
Kearny, NJ 07032-0607 USA  
Tel: 201.991.5330 Fax: 201.991.5870  
[www.LRultrasonics.com](http://www.LRultrasonics.com) [info@LRultrasonics.com](mailto:info@LRultrasonics.com)