



General Purpose Cleaner Non-Ammoniated

Removes blood, tissue, debris and other contaminants from instruments when used in an ultrasonic cleaning machine. Can also be used as a carrier bath for beakers and auxiliary pans. Concentrated. Available in gallons.

Product Code 228



General Purpose Cleaner Ammoniated

Extra strength for laboratory use. Ideal for instruments, molds, bridges and crowns. Brightens the appearance of both precious and non-precious alloys when used in an ultrasonic cleaning machine. Concentrated.

Available in gallons.

Product Code 226



Plaster and Stone Remover

Removes tough-to-clean plaster and stone from dentures, partials, flasks and articulators when used in an ultrasonic cleaning machine. Available in gallons.

Product Code 230



Tartar, Light Stain and Permanent Cement Remover

Removes tartar and light stains from dentures, partials and permanent cements when used in an ultrasonic cleaning machine. Can also be used to clean solder flux, oxide, investment, calcium hydroxide, some temporary cements and hard water scale. Available in gallons.

Product Code 232



Advanced Formula Temporary Cement Remover

Removes today's most popular temporary cements and replaces manual removal when used in an ultrasonic cleaning machine. Available in gallons.

Product Code 293



The Ultimate Infection Control Specialty Cleaning Solutions

When it comes to cleaning, the L&R reputation remains unsurpassed. For over 75 years, L&R has consistently provided leadership in the areas of cleaning innovations, safety, infection control and reliability and durability of product. L&R offers only the highest quality systems and premium solutions.

L&R solutions are specially formulated and rigorously tested for the ultimate cleaning capabilities, with each solution designed for a specific cleaning function. This is a time-honored tradition at L&R and is evidenced by the longevity of the products in the marketplace.

The practitioner has long associated L&R with the highest level of performance – performance recognized worldwide.

AutoClean

Keeps autoclaves operating at maximum performance

L&R's AutoClean removes hard water scale and lime deposits that accumulate in the chamber and helps maintain the efficiency of an autoclave. Available in quarts.

Product Code 239



PreZyme

Enzymatic cleaning formula

PreZyme Pre-Soak immediately begins removing proteinaceous soils, such as dried blood, mucous, debris and foreign matter and plays an important part in maintaining proper asepsis. Multi-use formula can also be used as an ultrasonic cleaning solution or to clean evacuation systems. Available in gallons.

Product Code 242



Barrier Milk

Protects your instruments during autoclaving

For use after ultrasonic cleaning and prior to autoclaving. A one-minute dip in Barrier Milk provides comprehensive protection from rust, heat discoloration and corrosion and keeps hinged instruments and box locks moving smoothly. Available in gallons.

Product Code 076



E-VAC

Non-foaming formula eliminates bacteria-breeding debris

Specially formulated to help dissolve organic debris, rope saliva and blood, the powerful non-foaming solution also cleans away calcium buildup in the lines. Traps, pumps, drains and lines are left minty fresh, without manual cleaning. Available in gallons.

Product Code 107



Call your representative today or log on to www.LRultrasonics.com for more information on our Ultrasonic Cleaning Solutions and entire product line.