



SweepZone[®] Ultrasonic Cleaning Systems

All SweepZone Ultrasonic Cleaning Systems operate with a special cleaning wave of +/- 2KHz energy that literally sweeps back and forth throughout the ultrasonic tank. The innovative L&R SweepZone technological system incorporates Pulse Width Modulation (PWM), which enables the machines to sense cleaning loads, solution levels and voltage fluctuations.

SweepZone Ultrasonic Cleaning Systems automatically adjust to different conditions in order to create more uniform cleaning power by eliminating hot spots and weak areas in the tank. This results in unsurpassed cleaning performance with every application.



SweepZone 200

User friendly digital touch pad controls make it easy and convenient to set any of the five preset timer settings from a range of six to sixty minutes. Durable stainless steel housing and tank are corrosion resistant and long lasting.

SweepZone® automatically adjusts for more uniform cleaning.



SweepZone 310



SweepZone 310 Recessed System
(Shown with hinged lid upgrade)



SweepZone 1000

Model*	Input Power Average watts 117 VAC 60 Hz	Amps	Overall Dimensions** Inches & Centimeters (L x W x H)	Shipping Weight Pounds & Kilograms	Tank Capacity Gallons & Liters	Tank Internal Dimensions Inches & Centimeters (L x W x D)
SweepZone 200	100	0.812	12.68 x 6.88 x 8.75 32.2 x 1.75 x 22.2	12 5.4	1.06 4.0	11.75 x 6 x 4 29.8 x 15.2 x 10.2
SweepZone 310	400	3.419	16.5 x 10 x 2 41.9 x 25.4 x 30.5	26 11.8	3.25 12.3	15.5 x 9 x 6 39.4 x 22.9 x 15.2
SweepZone 310R (Recessed)	400	3.418	16.25 x 10.13 x 6.25 41.3 x 25.7 x 15.9	26 11.8	3.25 12.3	15.5 x 9 x 6 39.4 x 22.9 x 15.2
SweepZone 360R (Recessed)	300	2.564	16.25 x 10.13 x 6.25 41.3 x 25.7 x 15.9	22 10.0	3.59 13.6	11.5 x 9.25 x 8 29.2 x 23.5 x 20.3
SweepZone 650R (Recessed)	500	4.274	16.25 x 10.13 x 6.25 41.3 x 25.7 x 15.9	38 17.2	6.50 24.6	19.75 x 11.75 x 6.5 50.2 x 29.8 x 16.5
SweepZone 1000	500	4.274	17 x 19 x 20 43.2 x 48.3 x 50.8	101 45.8	11.62 44.0	14 x 16 x 10 35.6 x 40.6 x 25.4

*All models include: Timer, Drain & Stainless Steel Cover

**Overall Dimensions of Recessed Systems are for Generator only

All systems operate at 43kHz - 45kHz sweep range with a sweep rate of either 3 or 83Hz. All systems are available with international voltages.

All L&R SweepZone Ultrasonic Cleaning Systems are backed by an unprecedented
2 ½ year warranty on parts and labor.

Visit our website to view our full line of accessories

L&R Manufacturing Company
577 Elm Street, PO Box 607
Kearny, NJ 07032-0607 USA
Tel: 201.991.5330 Fax: 201.991.5870
www.LRultrasonics.com info@LRultrasonics.com

