

# SweepZone® Recessed Ultrasonic Cleaning Systems



(Cabinet Not Included)

L&R's Recessed Systems with digital touch pad controls easily install into existing countertops and work areas to save valuable workspace.

These models provide an efficient and convenient method for instrument cleaning & processing. SweepZone recessed systems utilize simple inline connections that facilitate installation with most cabinet and Steri-Center Work Stations.

Product Code	Description	Input Power	Output	Generator Dimensions	Shipping Weight	Tank Capacity	Tank Internal Dimensions
		Average watts 117 VAC 60 Hz	Frequency	Inches Centimeters (L x W x H)	Pounds Kilograms	Quarts & Liters	Inches Centimeters (L x W x D)
AG 916	<b>SweepZone 310R</b> w/Remote Digital Timer, Drain, Mounting Rim	400	43 KHz +/-2KHz 83Hz	16.25 x 10.13 x 6.25 41.3 x 25.7 x 15.9	26 11.7	13 12.3	15.5 x 9 x 6 39.4 x 22.9 x 15.2
AG 915	<b>SweepZone 360R</b> w/Remote Digital Timer, Drain, Mounting Rim	300	43 KHz +/-2KHz 83Hz	16.25 x 10.13 x 6.25 41.3 x 25.7 x 15.9	22 9.9	14.5 13.6	11.5 x 9.25 x 8 29.2 x 23.5 x 20.3
AG 917	<b>SweepZone 650R</b> w/Remote Digital Timer, Drain, Mounting Rim	500	43 KHz +/-2KHz 83Hz	16.25 x 10.13 x 6.25 41.3 x 25.7 x 15.9	38 17.2	26 24.6	19.75 x 11.75 x 6.5 50.2 x 29.8 x 16.5



Product Code	Stainless Steel Solid Side Draining Baskets
13518	310, 310R (Full Size)
18202	310, 310R (1/2 Size)
18676	360, 360R (Full Size)
17954	650, 650R (Full Size)



Solid Side Basket (PC 13518)



310 Cassette Rack (PC 26665)

Product Code	Instrument Cassettes & Racks
19368	360R - Cassette Rack (Holds 5 Small Cassettes)
17582	Bur Box (Holds 32 Burs)
26665	310, 310R - Cassette Rack (Holds 3 Large Cassettes or 6 Small Cassettes)
26847	650, 650R - Cassette Rack (Holds 3 Large Cassettes or 6 Small Cassettes) 2 Racks Fit 1 Unit

Product Code	Accessory Kits
454	140 - Solid Side Basket, Plastic Positioning Cover, 2-600 ML Beakers Complete, Bur Tray
535	310, 310R - 1/2 Size Solid Side Basket, 1/2 Size Positioning Cover, 2-600 ML Beakers Complete, Bur Box, Bur Tray
26366	200 - 1 600 ML Beaker, 1 Plastic Beaker Cover, Full Size Stainless Steel Basket, 1/2 Size Positioning Cover, Bur Box, Bur Tray



**L&R Manufacturing Company**  
 577 Elm Street, P.O. Box 607  
 Kearny, NJ 07032 0607 USA  
 (P) 866.567.8811 / 201.991.5330 (F) 201.991.5870  
 www.LRultrasonics.com info@LRultrasonics.com