



UltraDose®

ENZYME PLUS


Ultrasonic Cleaning Solution

Welcome to the power of 3 with UltraDose® Enzyme Plus Ultrasonic Cleaning Solution. Designed to work with an Ultrasonic Cleaning System, you'll get broad spectrum cleaning due to our blend of three unique enzymes and surfactants that targets and breaks down proteins, carbohydrates, and fats. Enzyme Plus removes organic debris, residues, biofilm, inorganic soils, tissue, mucous, and blood from instruments, and general medical equipment. Plus, it's safe on stainless steel, metals, glass, plastics, rubber, and aluminum. Enzyme Plus delivers excellent cleaning results every cleaning cycle.

Unleash the cleaning power of enzymes

- Key part of infection prevention and instrument reprocessing
- Super concentrated, 1oz. makes one gallon of RTU solution
- Measured, easy-to-use pour dispensing system
- Ph balanced and free rinsing formula
- A blend of three unique enzymes that delivers protein-targeted cleaning
- Can also be used as a presoak

THE POWER OF 3



32 oz. Bottle
Pr. Code: UD038

L&R Manufacturing Company
577 Elm Street, P.O. Box 607
Kearny, NJ 07032-0607 USA
Tel: 201.991.5330 Fax: 201.991.5870

www.LRultrasonics.com info@LRultrasonics.com

because clean matters.



UltraDose[®]

ENZYME PLUS

Ultrasonic Cleaning Solution



It's not truly clean until it's UltraDose[®] clean. UltraDose Enzyme Plus attacks and breaks down protein chains and suspends them for fast and easy ultrasonic cleaning. The result is a thorough cleaning you can count on every time. Don't trust your instruments with just any cleaner: UltraDose Enzyme Plus attacks protein chains on a molecular level. It's the Enzyme Plus difference from L&R. Target problems and get a thorough cleaning—L&R provides the 360° cleaning solution that professionals rely on when cleaning matters most.

PROTEIN CHAINS


Protein chains including organic debris, blood, and mucous are formed on instruments.

BREAK DOWN

Enzymatic action attacks and breaks down tough protein chains.

SUSPEND AND REMOVE

Proteins are suspended in ultrasonic cleaning solution for easy removal.



because clean matters.

L&R Manufacturing Company
577 Elm Street, P.O. Box 607
Kearny, NJ 07032-0607 USA
Tel: 201.991.5330 Fax: 201.991.5870

www.LRultrasonics.com info@LRultrasonics.com

ISO 9001:2015

SOL.101.0518R