

UltraDose®

The UltraDose family of products gives you the ultimate in convenience and cost-effectiveness. These single-dose, super-concentrated powders and solutions require less storage and mixing and replace tablets and bulky packaging. The powerful new cleaning formulations of UltraDose products are easy to use, eliminate inaccurate measuring and dissolve quickly. Contact your dealer representative today for UltraDose products.



Dealer



www.LRultrasonics.com

phone: 201-991-5330 • fax: 201-991-5870

email: info@LRultrasonics.com

577 Elm Street, PO Box 607,

Kearny, New Jersey 07032-0607 USA

300-17-705-REV1-2.5M-STJ



Ultrasonic Cleaning Solution

An excellent cleaner of instruments, burs and orthodontic appliances, this solution removes blood, organic debris, buffing compounds, soils, oxides, rouges and tripoli in an ultrasonic cleaning machine. Each 1 ounce packet makes one gallon of cleaning solution. Each box makes 24 gallons of cleaning solution.

Product Code 012



General Purpose Cleaner

This powder cleans instruments, burs, general equipment, orthodontic appliances, gold and porcelain. Effectively removes blood, tissue, debris, buffing compounds, rouge, tripoli and oxides when used in an ultrasonic cleaning machine. One packet makes one gallon of cleaning solution. Each box makes 24 gallons of cleaning solution.

Product Code UD030

Tartar & Light Stain Remover

This powder cleans dentures, partials and bridges when used with a beaker in an ultrasonic cleaning machine. Biodegradable. Use one packet per cleaning. Each box contains 24 applications.

Product Code 028



Plaster & Stone Remover

This powder effectively removes heavy deposits of plaster and stone from dentures, partials, flasks, articulators and general instruments when used in a beaker or auxiliary pan in an ultrasonic cleaning machine. Use one packet per cleaning. Each box contains 12 applications.

Product Code UD031

UltraDose Germicidal

This Hospital Grade Ultrasonic Cleaner Concentrate cleans and disinfects inanimate, non-porous instrument and appliance surfaces, eliminating residual micro-organisms in one step. Performs as an ultrasonic cleaner, disinfectant, deodorant, germicidal, fungicidal, virucidal and tuberculocidal all-in-one. Each bottle makes 16 gallons of non-irritating, fragrance-free solution. (USA only)

Product Code UD036

