


## Products Rated Highly by Evaluators in CR Clinical Trials

### UltraDose WaveCheck



*L & R Mfg. Company*  
(sold through distributors)

**\$114/Box of 50**  
(~\$2/Indicator)

### Ultrasonic Cleaner Indicator Provides Visual Assurance of Adequate Instrument Debris Removal

Ultrasonic cleaners are helpful cleaning devices that should be tested for effectiveness, but frequently are not. This ultrasonic cleaner indicator was designed for testing/monitoring the effectiveness of ultrasonic cleaner performance. It is a ready-to-use monitor with red color on the surface of a tab to mimic blood and visually confirm debris removal and accuracy of device. They are helpful in assuring the effectiveness of ultrasonic cleaners.

#### Advantages:

- Easy to use in multiple locations of ultrasonic tank
- Results are quick, convenient, and easy to interpret
- Provides assurance that ultrasonic cleaner is working

#### Limitation:

- Some Evaluators desired an attachment/clip for the tab that is more adjustable.

**CR CONCLUSIONS:** 90% of 20 CR Evaluators stated they would incorporate UltraDose WaveCheck into their practice. 95% rated it excellent or good and worthy of trial by colleagues.